

























































# SURVEY OF CAPILLARY-PROFILES

## Capillary profiles (single)

	single, circular, thin-walled		single, circular, with hollow filament		theta-profile with filaments on inner wall		oval-capillary-profile		six-lobe profile, circular inside, without filament
	single, circular, standard wall		single, circular, with solid filament		theta cross-section-profile		2-hole-capillary-profile		three-lobe profile, circular inside, without filament
	single, circular, thick-walled		theta-profile, circular, with internal septum		single, square, standard wall thickness		4-hole-capillary-profile		triangular-profile with hole
	single, circular, very thick-walled		theta profile with filaments on internal septum		flat-capillary-profile		internal profile, circular outside, without filament		triangular-capillary-profile

## Capillary profiles (multiple)

	dual, different diameters, with solid filaments		triple bundle, same diameters, without filaments		quadruple bundle, same diameters, with filaments		sevenfold bundle, same diameters, without filaments		triple bundle, different diameters, with flanking rods
	dual, same diameters, without filaments		triple bundle, same diameters, with only one filament		quadruple bundle, same diameters, without filaments		sevenfold bundle, same diameters, with only 6 filaments		1+4er internal geometry, circular, without filaments
	dual, same diameters, with only one filament		triple bundle, same diameters, with only 2 filaments		quintuple bundle, same diameters, with filaments		sevenfold bundle, same diameters, with hollow filaments		multi-channel capillary-profile
	dual, same diameters, with filaments		triple parallel, same diameters, with filaments		quintuple bundle, same diameters, with only 4 filaments		3+4er geometry, bundled, with filaments		multi-channel capillary-profile (circular)
	dual, same diameters, with hollow filaments		triple parallel, same diameters, with only 2 filaments		quintuple bundle, same diameters, without filaments		quintuple bundle, same diameters, in x-array with hollow filaments		
	piggybag-profile		triple parallel, same diameters, without filaments		quadruple parallel, same diameters, without filaments		triple bundle, different diameters, with filaments outside		
	triple bundle, same diameters, with filaments		triple parallel, same diameters, with center rod		quadruple parallel, same diameters, with filaments		quintuple bundle, different diameters, in x-array with only 2 filaments		
	triple bundle, same diameters, with hollow filaments		triple, same diameters, with only one filament		sevenfold bundle, same diameters, with filaments		triple bundle, different diameters, with flanking capillaries, without filaments		